

Introduction to TMS

TMS Basics

Alvaro Pascual-Leone, MD, PhD

Senior Scientist, Hinda and Arthur Marcus Institute for Aging Research
Medical Director, Deanna and Sidney Wolk Center for Memory Health, Hebrew SeniorLife
Professor of Neurology, Harvard Medical School, Boston MA, USA



Marcus Institute
for Aging Research
Hebrew SeniorLife





HARVARD MEDICAL SCHOOL
AFFILIATE

TMS Basics


- Misnomer
- Stimulation to brain interaction
 - Brain structure
 - Brain state
- Pulse shape
 - Single pulse
 - Pairs of pulses
 - Trains of repetitive pulses
- Variability


Electro-Magnetic Induction

“I think I got hold of a good thing”




M. Faraday
29 August 1831





Anthony Barker
1984



Origin of Therapeutic TMS – Repetitive TMS

1984




Anthony Barker
Single Pulse TMS

1987




Cadwell
Repetitive TMS
(rTMS)

1996



Pascual-Leone, et. al.
George, et. al.
rTMS for depression




Pascual-Leone et al. The Lancet 1996

2007



Neuronetics Phase III
trial of rTMS for
Medication-resistant
depression



FDA
approval

2008

NHIC
Medicare
Approval
(MA, NH, VT and
RI)

2012

Nexstim

FDA approved for cortical mapping (motor & speech)




Umer
Nagib

Shahid
Bashir


Alex
Golby

HAND
(APB, ADM)

Foot
(TA)



Tumor




Tumor

Craniotomy

Foot
(TA)

FDA Cleared Neuronavigation



Robot-guided TMS

- Target different network nodes in specific order and timing
- Adaptive close loop system



FDA Cleared Devices for
 Medication-Resistant Depression
 Adolescent depression
 Migraine
 OCD
 Substance Abuse D/O
 Anxiety D/O

Growing list of off-label applications

PTSD
 Bipolar disorder
 Autism
 Pain
 Parkinson's Disease
 Focal Epilepsy
 Impulse control D/O
 Agitation
 MCI & Dementia



magstim

eNeura[®]



Brainsway

MagVenture



Nexstim



Real Clinical Impact !

TMS in Medication-Resistant Depression

- >600 systems in clinical use in the US
- 250 days/year & 5 patients/day
= *750,000 treatments per year*
- approx. 25 sessions/Rx/patient
= *30,000 patients/year*
- 30% remission
= *9,000 patients in remission/year*
- *25 patients in remission/day*

Covered by Medicare
& most health insurance plans in the US

Covered by health insurance in

- Canada
- Australia
- New Zealand
- Japan
- UK


Real Clinical Impact !

TMS in Medication-Resistant Depression

- After a treatment course (of up to 6 weeks) benefit lasts on average 5 months
- In case of relapse, response to new treatment course is at least as good as initial response in >90% of cases [Kelly et al. J Neuropsych Clin Neurosci 2017]
- Maintenance is possible

Helping Patients
However.....only the beginning

Realize the Promise of Neuromodulation



Noninvasive Neuromodulation does not represent a treatment for an illness, but offers **tool** that allows modulation of the neural substrate of **symptoms and disabilities** caused by brain illnesses or dysfunctions



Helping Patients

However.....only the beginning

Realize the Promise of Neuromodulation

PRECISION

- Competencies & Training
- Safety
- Spatio-temporal signatures of brain-related disability
- Target specific symptoms/disabilities
- Leverage individual differences

IMPORTANCE OF MODELLING • NEW TECHNOLOGIES • SCALABILITY

- Home
- New Technologies
- Targeting specific symptoms/disabilities





ELSEVIER



Guidelines

Safety, ethical considerations, and application guidelines for the use of transcranial magnetic stimulation in clinical practice and research[☆]

Simone Rossi^{a,*}, Mark Hallett^b, Paolo M. Rossini^{c,d}, Alvaro Pascual-Leone^e and The Safety of TMS Consensus Group¹



ELSEVIER



Review

Safety and recommendations for TMS use in healthy subjects and patient populations, with updates on training, ethical and regulatory issues: Expert Guidelines

Simone Rossi^{a,*}, Andrea Antal^{b,c}, Sven Bestmann^d, Marom Bikson^e, Carmen Brewer^f, Jürgen Brockmüller^g, Linda L. Carpenter^h, Massimo Cincottaⁱ, Robert Chen^j, Jeff D. Daskalakis^k, Vincenzo Di Lazzaro^l, Michael D. Fox^{m,n,o}, Mark S. George^p, Donald Gilbert^q, Vasilios K. Kimiskidis^r, Giacomo Koch^s, Risto J. Ilmoniemi^t, Jean Pascal Lefaucheur^{u,v}, Letizia Leocani^w, Sarah H. Lisanby^{x,y,2}, Carlo Miniussi^z, Frank Padberg^{aa}, Alvaro Pascual-Leone^{ab,ac,ad}, Walter Paulus^b, Angel V. Peterchev^{ae}, Angelo Quartarone^{af}, Alexander Rotenberg^{ag}, John Rothwell^d, Paolo M. Rossini^{ah}, Emiliano Santarnecchi^m, Mouhsin M. Shafi^m, Hartwig R. Siebner^{ai,aj,ak}, Yoshikatzu Ugawa^{al}, Eric M. Wassermann^{am,2}, Abraham Zangen^{an}, Ulf Ziemann^{ao}, Mark Hallett^{ap,2,*}

The basis of this article began with a Consensus Statement from the IFCN Workshop on “Present, Future of TMS: Safety, Ethical Guidelines”, Siena, October 17–20, 2018, updating through April 2020¹



ELSEVIER



Review

Training in the practice of noninvasive brain stimulation: Recommendations from an IFCN committee


Peter J. Fried^{a,1}, Emiliano Santarnecchi^{a,1}, Andrea Antal^b, David Bartres-Faz^c, Sven Bestmann^d, Linda L. Carpenter^e, Pablo Celnik^f, Dylan Edwards^{g,h}, Faranak Farzanⁱ, Shirley Fecteau^j, Mark S. George^{k,l}, Bin He^m, Yun-Hee Kimⁿ, Letizia Leocani^o, Sarah H. Lisanby^p, Colleen Loo^q, Bruce Luber^r, Michael A. Nitsche^s, Walter Paulus^b, Simone Rossi^t, Paolo M. Rossini^u, John Rothwell^v, Alexander T. Sack^w, Gregor Thut^x, Yoshikazu Ugawa^y, Ulf Ziemann^z, Mark Hallett^{aa}, Alvaro Pascual-Leone^{ab,ac,*}




Realize the Promise of Neuromodulation

PRECISION


Spatial precision




Temporal precision



Contextual precision






Spatial precision




1. Know where to target
 - MRI-guided TMS
2. Keep target consistent
 - Robot-assisted TMS
3. Make target smaller
 - Micro TMS



NiBS modulates activity in brain networks & the effects depend on connectivity



Functional: TMS impact (% Δ ^{14}C 2-DG)




- A19
- A18
- A17
- SVA
- CVA
- pCG
- PUL
- LPI
- ▲ SGS
- ▲ SO

Anatomy: Projection Weight (grains/ μm^2)

Impact on specific brain networks




Davide Momi






Network Engagement

STIMULATION over DMN

STIMULATION over DAN





The effect of coil orientation





Giorgio Bonmassar

μ TMS










Mike Fox

Brain Target and Need to Individualize

rs-fcMRI connectivity




Clinical Efficacy vs.
Subgenual Ncl Anticorrelation




Temporal precision


Waveform optimization



Patterned stimulation



Closed-loop stimulation






1. Individualize stimulation parameters
 - Measure neurophysiologic effect
2. Design improved stimulators
3. Stimulate at right time
 - Closed loop stimulation




Antoni Valero-Cabré


rTMS modulates activity in targeted brain beyond the duration of the stimulation & depending on stimulation frequency





Sham TMS



1 Hz TMS



20 Hz TMS




Variability of Physiologic Effects



Fumiko
(Maeda) Hoeft

Δ % in MEP area
pre/post rTMS




There is no
'inhibitory' or
'excitatory' rTMS


Variability of Physiologic Effects





Mark Eldaief




Low frequency rTMS






High frequency rTMS



Need to Measure!
Neurophysiologic monitoring:
fMRI - EEG - etc
Define dose
Enable close-loop




Neurophysiologic Biomarkers Predict Clinical Response




Contextual precision

Online stimulation



Biological rhythm



Combinatorial therapy



1. Stimulate at right time

- Closed loop stimulation

2. Leverage State-dependency effects





Zaira
Cattaneo



Juha
Silvanto

State Dependency For Specific Neural Populations

TMS-adaptation paradigm and motion direction discrimination in visual area V5/MT



TMS disrupts non-adapted but improves adapted direction discrimination


State-Dependency of Transcranial Magnetic Stimulation

Juha Silvanto · Alvaro Pascual-Leone

Concurrent TMS stimulation with Cognitive Training in Dementia




FDA-Approval for OCD
Provocation task + TMS



Brainsway

- How create an 'optimal' state ?
- When to couple brain stimulation with state modification ?



Factors Contributing to Individual Variability

- Genetics
 - BDNF met/val
- Age
 - Age-related GABA effect
- Head shape
 - Asian vs Caucasian
- IQ, Cognitive Capacity
- Brain structure
 - Atrophy
 - Connectivity
- Circadian rhythms
 - Time of day
- Hormonal state
 - Menstrual cycle
- Metabolic state
 - Low K⁺ (diarrhea)
- Nutrition considerations
 - Coffee
- Medications
 - Valproate
- State / Mindset

Leverage Individual Differences
instead of trying to Control them?

TMS Basics

- Misnomer
 - Stimulation to brain interaction
 - Brain structure
 - Brain state
 - Pulse shape
 - Single pulse
 - Pairs of pulses
 - Trains of repetitive pulses
 - Variability
- Electrodeless Electric Stimulation via Electromagnetic Induction
 - There is no FOCAL TMS – CIRCUIT THERAPEUTICS
 - Not a treatment for a disease – but for symptoms/disabilities
 - Cortical target as ‘window’ – leverage brain network knowledge
 - There is no INHIBITORY or EXCITATORY TMS
 - There are huge interindividual differences
 - Need to measure – relate behavioral to physiologic effects
 - Effects depend on brain state – STATE DEPENDENCY is critical
 - Use state dependency for greater effect specificity
- REMEMBER:
- Details matter (coil size and shape, pulse shape, current direction, etc.)
 - More is not necessarily better – NEUROPHYSIOLOGIC PRECISION
 - Test-rest reliability is poor and should be assessed
 - Sham stimulation is tricky – ACTIVE controls